

Cellacare® Genucast Modular

Knee splint for posttraumatic and postoperative immobilisation of the joint in a functionally bent position of 15° dorsal in cases of inflammation and injuries of various types.

Removable aluminium lateral and dorsal stays secure and stabilise the injured knee. Cellacare Genucast Modular features hook & loop closures which are quick and easy to use. This orthosis consists of three elements which makes it possible to fit any individual leg circumference. Open-cell foam padding is lined with skin-friendly cotton material providing high wearer comfort.

Cellacare® Genucast Modular

Foam padding: polyurethane/polyester; white inside lining: cotton, blue outside cover: cotton, edging: cotton, hook & loop tapes and eyelets: polyamide, longitudinal tapes: polypropylene, sewing thread: cotton, polyester, aluminium stays.

Washable up to 30°C, remove stays before washing, do not iron, do not tumble-dry, not suitable for dry cleaning, colour: navy/grey.



Cellacare® Genucast Modular

individually in PE bag

Size	for Body Size	Height	REF	Shipping Units
1	< 1.60 m	43 cm	22 472	1 pc. / 14 pcs.
2	< 1.60 bis 1.80 m	53 cm	22 473	1 pc. / 10 pcs.
3	> 1.80 m	63 cm	22 474	1 pc. / 10 pcs.