

Cellacare® Gilchrist

Cellacare® Gilchrist Easy

Ready-to-wear supports on the principle of the Gilchrist or De-sault bandage for the posttraumatic or postoperative immobilisation of the shoulder and upper arm. Also for severe distortions or fractures in the region of the shoulder and upper arm, rupture of the acromioclavicular joint, rotator injuries, shoulder luxation and impingement syndrome.

The ready-to-wear supports can be put on with little manipulation and in a very short time, causing neither pain nor pressure, because the injured shoulder remains free.

Therapeutic exercises and wound dressing can be undertaken without removing the support. Daily hygiene will scarcely be affected.

Flat, padded straps made from skin-friendly material.

The practical handling of these supports provides a high level of wearer comfort.

Cellacare® Gilchrist

The hook & loop fasteners are fitted in such a way that the support can be used on the left or the right. 8 cm wide waistband: foam padding polyurethane, outer cover: viscose, hook & loop tape: polyamide, velour tapes: polyamide quilted seams: PES/cotton.

Washable: up to 30°C, does not require ironing, not suitable for dry-cleaning, colour: ecru.

Cellacare® Gilchrist Easy

Support for the immobilisation of the shoulder joint in a functional position; reusable, ready-to-wear, 100% cotton, skin-friendly, washable up to 95°C, colour: ecru.



Cellacare® Gilchrist



Cellacare® Gilchrist Easy

Cellacare® Gilchrist

individually packed in a cardboard box

Size	Length	Chest Size	REF	Shipping Units
S	190 cm	up to 60 cm	22 401	1 pc. / 44 pcs.
M	230 cm	60- 80 cm	22 402	1 pc. / 44 pcs.
L	280 cm	80-100 cm	22 403	1 pc. / 44 pcs.
XL	320 cm	100-120 cm	22 404	1 pc. / 44 pcs.

Cellacare® Gilchrist Easy

individually packed in a cardboard box

left	1	approx. 80-100 cm	22 415	1 pc. / 5 pcs.
	2	approx. 100-120 cm	22 416	1 pc. / 5 pcs.
right	1	approx. 80-100 cm	22 417	1 pc. / 5 pcs.
	2	approx. 100-120 cm	22 418	1 pc. / 5 pcs.