

# Cellacare® Lumbal F

# Cellacare® Lumbal M

Anatomically shaped back support, available in two versions, for the male and female anatomy. With a perfect fit, individually adjustable, available in various cuts and 5 sizes each. Rigid stays in the back and stabilising stays in the stomach region ensure a safe stabilising performance. The individual adjustment of the lumbar stability and the pressure of the insert is achieved with the transverse straps of the back support. The innovative, ergonomically shaped WaveTec Insert which is safely attached to the support with a hook & loop strip, improves the neuromuscular control in the region of the back

muscles. The design with ergonomically rounded edges and a specific, breathable material composition with a high cotton content offers excellent wearer comfort.

Cellacare Lumbal relieves the intervertebral disks when lifting and carrying weights, stabilises the lumbar spine and corrects false posture. Thereby, Cellacare Lumbal reduces pain, assists treatment and has a preventive function by changing patterns of behaviour.

Indications: for the treatment and prevention of acute and chronic lumbalgia, lumbar ischialgia and lumbar syndrome,



protrusions, prolapse, facet syndrome and spondylolysis, chronic false posture and irritations.

### Cellacare® Lumbal

Support: cotton, polyamide, polyester, elastodiene, polypropylene. Insert: polyurethane foam.

Washable up to 30°C, do not iron, do not tumble-dry, not suitable for dry-cleaning, colour: grey/dark grey.



Cellacare® Lumbal F



Cellacare® Lumbal M

### Cellacare® Lumbal F

Back support for women

Size	Waist Circumference	REF	Shipping Units
1	60 - 75 cm	13 368	1 pc. / 12 pcs.
2	75 - 95 cm	13 369	1 pc. / 12 pcs.
3	95 - 115 cm	13 370	1 pc. / 12 pcs.
4	115 - 135 cm	13 371	1 pc. / 12 pcs.
5	> 135 cm	13 372	1 pc. / 12 pcs.

### Cellacare® Lumbal M

Back support for men

1	60 - 75 cm	13 373	1 pc. / 12 pcs.
2	75 - 95 cm	13 374	1 pc. / 12 pcs.
3	95 - 115 cm	13 375	1 pc. / 12 pcs.
4	115 - 135 cm	13 376	1 pc. / 12 pcs.
5	> 135 cm	13 377	1 pc. / 12 pcs.