

Cellacare® Malleo

Active ankle support, for compression and support of the ankle joint, with a moulded heel for perfect anatomical fit; for the functional treatment of post-traumatic and postoperative conditions, distortion, contusion, tendopathy, ligament instability, intra-articular effusion, and swelling caused by active arthrosis or arthritis, as well as for injury prevention.

Firm on the outside and soft on the inside - due to the specific material composition and three-dimensional knitting.

The moulded heel, additional pressure-relief zones on the heel and back of the foot, and the inserts sewn firmly onto the inner and outer ankle sections provide an improved fit, and thus a better distribution of pressure. The special

massaging effect, which increases mobility and stimulates the metabolism, is also efficient in reducing swelling and pain.

Cellacare Malleo has a stabilising and compressing effect on the ankle joint, thus improving and re-adjusting the alignment of the joint. The inserts ensure a physiologically appropriate distribution of pressure. The support is suitable for both the right and left foot and can be worn with all common walking shoes or sneakers.



Cellacare® Malleo

Three-dimensional knitting,
colour: blue/black,
cotton/nylon fabric,
inserts: silicone,
insert pockets: cotton.



Soft cotton inside

Washable at temperatures of up to 30 °C,
do not tumble-dry,
do not iron,
not suitable for dry cleaning.

Cellacare® Malleo

individually packed in a cardboard box

Size	Ankle Circumference	REF	Shipping Units
1	17 - 19 cm	16 856	1 pc. / 20 pcs.
2	19 - 21 cm	16 857	1 pc. / 20 pcs.
3	21 - 23 cm	16 858	1 pc. / 20 pcs.
4	23 - 26 cm	16 859	1 pc. / 20 pcs.
5	26 - 29 cm	16 860	1 pc. / 20 pcs.
6	29 - 32 cm	16 861	1 pc. / 20 pcs.