

Cellacare® Malleo Akut

Orthosis for the stabilisation and protection of the ankle joint, according to the accepted classical U-splint mechanism of action; Main indications: fresh, fibular ligament ruptures and severe distortions of the lateral ligaments.

Cellacare Malleo Akut provides both medial and lateral enclosure of the ankle joint and prevents any twisting and thereby rupture or strain of the lateral

ligament apparatus. High wearer comfort through anatomically shaped shells with wide strap closures and soft stirrup. The inside padding is comfortable and skin-friendly, with an anti-bacterial cover.

Fitting instructions on the product allow easy application. The orthosis is comfortable to wear inside shoes. The stirrup can be individually adjusted with the product in place.



Cellacare® Malleo Akut

Side shells: polyethylene; Inside foam padding: polyester, polyurethane; Padding – textile cover: polyester microfibres, treated against bacteria/ polyamide; Stirrup: polyamide, polyester; Accessories: straps, buckles: polyamide; Fastening points: polypropylene.

Washable up to 30°C, do not iron, do not tumble-dry, not suitable for dry-cleaning, colour: blue/black.



Cellacare® Malleo Akut

individually in a cardboard box

Size	for Body Size	REF	Shipping Units
right			
1	< 1.60 m	22 466	1 pc. / 21 pcs.
2	> 1.60 m	22 468	1 pc. / 18 pcs.
left			
1	< 1.60 m	22 467	1 pc. / 21 pcs.
2	> 1.60 m	22 469	1 pc. / 18 pcs.