

Cellacare® Materna

The back support for expectant mothers, with flexible stays in the back and lateral pulls stabilises the abdomen, back, and pelvic region; assists circulation and activates muscles.

The relief of the back and pelvic floor and the simultaneous support of the abdomen reduce pain.

The wide elastic strap at the back supports the lumbar spine;

flexible adjustable straps are anatomically shaped and allow an individual adjustment of width. They are closed below the abdomen and provide comfort and support.

Indications: back complaints and signs of overexertion during and after pregnancy, for the support in cases of symphyseal disruption and instabilities after childbirth.



Cellacare® Materna

Basic textile: polyester, lycra; Side straps: nylon, cotton, PU; Elastic: polyester, lycra; Hook/loop fastener: nylon; Cover of stays: polyester; Thread: polyester/nylon; Edging: cotton, polyester. Colour: white.

Close all hook & loop closures before washing. Machine washable at 30°C. Do not dry in direct heat (heater, sun) or dryer.

Cellacare® Materna

individually packed in a cardboard box

Size	Hip Circumference	REF	Shipping Units
1	< 91 cm	17 008	1 pc. / 20 pcs.
2	91 - 102 cm	17 009	1 pc. / 20 pcs.
3	102 - 117 cm	17 010	1 pc. / 20 pcs.
4	> 117 cm	17 011	1 pc. / 20 pcs.