

Cellacare® Rhizocast

Thumb-wrist orthosis with a breathable 3D knit material as well as two volar aluminium stays and a thermoplastic insert in the area of the thumb for individual stabilisation.

Cellacare Rhizocast immobilizes the thumb in a functional position of abduction and opposition and, with its integrated reinforcements,

protects against external physical influences. Movement of the metacarpophalangeal joints (MCP joints) and the gripping function between thumb and finger tips are not inhibited.

The thumb-wrist orthosis is equipped with an application aid and closing straps with alternately placed buckles for simple, tangle-free fitting.

Cellacare Rhizocast is used for immobilization in cases of posttraumatic or postoperative irritations in the wrist area and the basal and saddle joints of the thumb, in cases of rhizarthrosis, ulnar ligament lesion (skier's thumb), STT arthrosis, tendovaginitis, carpal-canal syndrome, and distortion.



Cellacare® Rhizocast

Materials: polyester, polyurethane, cotton, elasthane, polyamide, polyvinyl chloride, polyethylene, aluminium, steel.

Colour: black-blue

Washable up to 30°C, do not tumble-dry, iron or dry-clean, no chlorine bleach.

Cellacare® Rhizocast

individually packed in a cardboard box

Size	Wrist Circumference	REF	Shipping Units
left			
1	< 17 cm	19 801	1 pc. / 20 pcs.
2	17 - 24 cm	19 802	1 pc. / 20 pcs.
3	> 24 cm	19 803	1 pc. / 20 pcs.
right			
1	< 17 cm	19 804	1 pc. / 20 pcs.
2	17 - 24 cm	19 805	1 pc. / 20 pcs.
3	> 24 cm	19 806	1 pc. / 20 pcs.