



Cellacare® Thorax F

Cellacare® Thorax M

Ready-to-wear supports for the conservative treatment of rib fractures and contusions of the thorax; also for postoperative relief after surgical interventions in the thoracic and abdominal region.

The two versions, M and F, take into account differences between the male and female anatomy.

The wide hook & loop closing guarantees correct fit and even distribution of pressure. The strong belt material with high transverse stability and 100% longitudinal stretch ensures even surface pressure distribution without ripples forming or the edges rolling over.

The padded ventral insert provides extra softness in the wound or injury area.



Cellacare® Thorax F

Cellacare thoracic supports can also be fitted satisfactorily on a supine patient.



Cellacare® Thorax M

Cellacare® Thorax

Foam padding: polyurethane/polyester, covers: inside cotton, outside polyester, elastic band: polyamide, cotton, elastodien, hook & loop: polyamide, edging: cotton, sewing thread: PES/cotton.

Washable up to 30°C, do not iron, do not tumble-dry, not suitable for dry-cleaning, colour: white.

Cellacare® Thorax F

individually packed in a cardboard box

Size	Thorax Circumf.	Height	REF	Shipping Units
S	70 to 100 cm	16 cm	22 470	1 pc. / 60 pcs.
L	100 to 130 cm	16 cm	22 471	1 pc. / 60 pcs.

Cellacare® Thorax M

individually packed in a cardboard box

S	70 to 100 cm	24 cm	22 475	1 pc. / 24 pcs.
L	100 to 130 cm	24 cm	22 476	1 pc. / 24 pcs.