

Ankle Supports

epX® Ankle Basic

Ankle support with open heel, anatomically sewn.

epX Ankle Basic supports and eases the upper ankle joint after injuries, e.g. distortions and contusions as well as in cases of acute and chronic pain caused by strain or overexertion. Also for degenerative joint diseases.



epX® Ankle Basic

epX® Ankle Dynamic

with application aid. Ankle support with stabilising strap, anatomically sewn. J-shaped inserts in the ankle region accelerate the breakdown of swellings and haematoma. The sewn-on strap is passed around the ankle in a figure of eight and secured with hook & loop closure.

epX Ankle Dynamic supports and eases the upper ankle joint after distortions, contusions,

capsular ligament lesions; also postoperatively (e.g. following suture of fibular ligament rupture) as well as in rehabilitation; also in cases of chronic ligament instability. Inserts are sewn into their pockets loosely, so that they will move automatically into the correct position around the inner and outer malleolus.

For details on materials and care see page 121.



epX® Ankle Dynamic

epX® Ankle Basic individually packed in a transparent plastic box				
Size		Ankle Circumference	REF	Shipping Units
S		17.5 - 20.5 cm	22 710	1 pc. / 24 pcs.
M		20.5 - 23.0 cm	22 711	1 pc. / 24 pcs.
L		23.0 - 25.5 cm	22 712	1 pc. / 24 pcs.
XL		25.5 - 28.0 cm	22 713	1 pc. / 24 pcs.
XXL		28.0 - 33.0 cm	22 714	1 pc. / 24 pcs.
epX® Ankle Dynamic individually packed in a transparent plastic box				
Right ankle version	S	17.5 - 20.5 cm	22 720	1 pc. / 24 pcs.
	M	20.5 - 23.0 cm	22 721	1 pc. / 24 pcs.
	L	23.0 - 25.5 cm	22 722	1 pc. / 24 pcs.
	XL	25.5 - 28.0 cm	22 723	1 pc. / 24 pcs.
	XXL	28.0 - 33.0 cm	22 724	1 pc. / 24 pcs.
Left ankle version	S	17.5 - 20.5 cm	22 725	1 pc. / 24 pcs.
	M	20.5 - 23.0 cm	22 726	1 pc. / 24 pcs.
	L	23.0 - 25.5 cm	22 727	1 pc. / 24 pcs.
	XL	25.5 - 28.0 cm	22 728	1 pc. / 24 pcs.
	XXL	28.0 - 33.0 cm	22 729	1 pc. / 24 pcs.