



# Back Support

## epX® Back Basic

Back support with elastic double straps and sewn-in stay above the spinal column for stabilisation. A large pocket is integrated in the dorsal area for an insert, which is also suitable for the use of thermoplastic inserts.

epX Back Basic supports the back during treatment and pre-

vention of acute and chronic lumbalgia or lumbal ischialgia. Also in cases of chronic false posture and irritations.

### Accessories for epX® Back Basic

To optimise the medical effect, inserts can be placed in the purpose-made pocket which is provided on the back support.

## Soft Insert

Conformable and soft, with a heating effect, made of coated foam material, with sewn-in insert; black, individually packed in a PE bag.



## Thermoplastic Insert

Made of thermoplastic material, can be heated in hot water and individually shaped to influence the pelvic tilt. (molding temperature minimum 50°C); white, individually packed in a PE bag.

For details on materials and care see page 121.

epX® Back Basic			
Size	Hip/Waist Circumference	REF	Shipping Units
S	86.0 - 96.5 cm	22 670	1 pc. / 12 pcs.
M	96.5 - 106.5 cm	22 671	1 pc. / 12 pcs.
L	106.5 - 117.0 cm	22 672	1 pc. / 12 pcs.
XL	117.0 - 127.0 cm	22 673	1 pc. / 12 pcs.
XXL	127.0 - 137.0 cm	22 674	1 pc. / 12 pcs.
Soft Insert			
universal		22 681	1 pc. / 30 pcs.
Thermoplastic Insert			
universal		22 680	1 pc. / 20 pcs.