

Elbow Supports

epX® Elbow Basic

Elbow support, sleeve type, with integrated insert. Its circular strap consists of an elastic and an inelastic part which is pulled through a buckle and simply attached with a hook & loop tape.

epX Elbow Basic supports the treatment and prevention of epicondylitis humeri radialis (tennis elbow) and ulnaris (golfer's elbow). The elbow support lessens the pain in cases of strain and overexertion and reduces the mechanical load at the base of the tendon.



epX® Elbow Basic

epX® Elbow Dynamic

Elbow support with insert and a circular strap over the insert, anatomically sewn.

epX Elbow Dynamic is indicated for the treatment and prevention of epicondylitis humeri radialis

(tennis elbow) and ulnaris (golfer's elbow). Also in cases of distortions and contusions, myotendopathies, as well as acute and chronic pain and signs of overexertion.

For details on materials and care see page 121.



epX® Elbow Dynamic

epX® Elbow Basic

individually packed in a transparent plastic box

Size	Forearm Circumference	REF	Shipping Units
S	22.5 - 25.5 cm	22 690	1 pc. / 60 pcs.
M	25.5 - 28.0 cm	22 691	1 pc. / 60 pcs.
L	28.0 - 30.5 cm	22 692	1 pc. / 60 pcs.
XL	30.5 - 35.5 cm	22 693	1 pc. / 60 pcs.
XXL	35.5 - 40.5 cm	22 694	1 pc. / 60 pcs.

epX® Elbow Dynamic

individually packed in a transparent plastic box

S	20 - 24 cm	22 700	1 pc. / 24 pcs.
M	24 - 26 cm	22 701	1 pc. / 24 pcs.
L	26 - 29 cm	22 702	1 pc. / 24 pcs.
XL	29 - 34 cm	22 703	1 pc. / 24 pcs.
XXL	34 - 40 cm	22 704	1 pc. / 24 pcs.