



Knee Supports

epX® Knee Basic

Knee support made of textile epX material with integrated lateral stays, anatomically sewn into a lightly shaped position with diagonal stretch in the thigh area; closed over the patella.

epX Knee Basic is indicated for the treatment and prevention of pain and overexertion of the knee joint. Also in cases of chronic pain or degenerative joint diseases.



epX® Knee Basic

This product is indicated for contusions, ligament lesions and chronic instabilities; also in the rehabilitation and prevention of patella-femoral pain syndromes, in cases of chondropathy patellae, myo- and tendopathies and overexertion. Also for the treatment of joint contusions and swelling caused by arthrosis and arthritis as well as for posttraumatic or postoperative conditions.

For details on materials and care see page 121.

epX® Knee Basic			
individually packed in a transparent plastic box			
Size	Knee Circumference	REF	Shipping Units
S	33.0 - 35.5 cm	22 600	1 pc. / 20 pcs.
M	35.5 - 37.0 cm	22 601	1 pc. / 20 pcs.
M/L	37.0 - 38.0 cm	22 602	1 pc. / 20 pcs.
L	38.0 - 40.5 cm	22 603	1 pc. / 20 pcs.
XL	40.5 - 45.5 cm	22 604	1 pc. / 20 pcs.
XXL	45.5 - 51.0 cm	22 605	1 pc. / 20 pcs.

Knee Supports

epX® Knee Dynamic

Knee support with locating straps, epX Knee Dynamic protects the damaged knee joint dynamically; lateral reinforcements and the two locating straps provide additional stabilization to the knee.

The ring-shaped insert can be cut to size and fitted into the open pocket of the support, thereby padding the patella.

epX Knee Dynamic is indicated for the treatment of contusions, ligament lesions and chronic instabilities (e.g. jumper's knee). Also in the rehabilitation and prevention of patella-femoral pain syndromes, in cases of chondropathy patellae, myo- and tendopathies and overexertion. Also for the treatment of joint contusions and swellings caused by arthritis and arthritis as well as for posttraumatic or postoperative conditions.



epX® Knee Dynamic



epX® Knee J Patella

epX® Knee J Patella

Knee support with integrated J-shaped insert. It has a centralising effect on the patella in cases of poor positioning, either inherent or due to injury. The lateral slipping of the patella during movement which causes the discomfort is effectively prevented through the J-shaped insert and the asymmetrically guided straps, as well as two lateral stabilising stays on each side. Five sizes for the right or left knee.

epX Knee J Patella is indicated for patella syndromes, chondropathy patellae, after patella luxation and for medialisation in cases of patellar deformation, also postoperatively e.g. after lateral release.

For details on materials and care see page 121.

epX® Knee Dynamic				
individually packed in a transparent plastic box				
Size		Knee Circumference	REF	Shipping Units
S		33.0 - 35.5 cm	22 620	1 pc. / 20 pcs.
M		35.5 - 37.0 cm	22 621	1 pc. / 20 pcs.
M/L		37.0 - 38.0 cm	22 622	1 pc. / 20 pcs.
L		38.0 - 40.5 cm	22 623	1 pc. / 20 pcs.
XL		40.5 - 45.5 cm	22 624	1 pc. / 20 pcs.
XXL		45.5 - 51.0 cm	22 625	1 pc. / 20 pcs.
epX® Knee J Patella				
individually packed in a transparent plastic box				
Right knee version	M	35.5 - 37.0 cm	22 630	1 pc. / 20 pcs.
	M/L	37.0 - 38.0 cm	22 631	1 pc. / 20 pcs.
	L	38.0 - 40.5 cm	22 632	1 pc. / 20 pcs.
	XL	40.5 - 45.5 cm	22 633	1 pc. / 20 pcs.
	XXL	45.5 - 51.0 cm	22 634	1 pc. / 20 pcs.
Left knee version	M	35.5 - 37.0 cm	22 635	1 pc. / 20 pcs.
	M/L	37.0 - 38.0 cm	22 636	1 pc. / 20 pcs.
	L	38.0 - 40.5 cm	22 637	1 pc. / 20 pcs.
	XL	40.5 - 45.5 cm	22 638	1 pc. / 20 pcs.
	XXL	45.5 - 51.0 cm	22 639	1 pc. / 20 pcs.