



Textile Orthosis

epX® Knee Control Plus

Textile orthosis with adjustable hinges. epX Knee Control Plus is provided with two circular straps on the thigh and lower leg and two straps immediately above and below the patella,

alternately wrapped in right and left direction, as well as lateral anatomically preshaped metal splints with hinged joints. For the protection of damaged ligaments, the range of movement during flexion and extension of the knee, which is guided by the

lateral splints, can be limited as required by adjusting the hinges in 10° and 20° steps.

A ring-shaped insert is integrated into the patella area.

epX Knee Control Plus guides the knee's mobility in a physiologically correct manner and provides additional security to unstable knee joints; for acute and chronic lateral ligament lesions, also postoperatively to limit knee mobility, as well as for the prophylaxis of new injuries during a rehabilitation phase.

For details on materials and care see page 121.



epX® Knee Control Plus

individually packed in a transparent plastic box

Size	Knee Circumference	REF	Shipping Units
S	33.0 - 35.5 cm	22 775	1 pc. / 12 pcs.
M	35.5 - 37.0 cm	22 776	1 pc. / 12 pcs.
M/L	37.0 - 38.0 cm	22 777	1 pc. / 12 pcs.
L	38.0 - 40.5 cm	22 778	1 pc. / 12 pcs.
XL	40.5 - 45.5 cm	22 779	1 pc. / 12 pcs.
XXL	45.5 - 51.0 cm	22 780	1 pc. / 12 pcs.