

# Wrist Supports

## epX® Wrist Basic

Wrist support, wrap-around type, with practical hook & loop closure. The support has a thumb opening and a soft, comfortable thumb strap.

epX Wrist Basic supports the treatment after distortions and contusions of the wrist and the metacarpal area, as well as in

cases of pain and overexertion of the distal forearm.

## epX® Wrist Dynamic

Wrist support with stabilising stays for dorsal and volar stabilisation, and with a circular supporting strap above the wrist with hook & loop closure; anatomically sewn, with a thumb opening and a soft thumb strap.



epX® Wrist Dynamic



epX® Wrist Basic

epX Wrist Dynamic supports and relieves the wrist and the metacarpal area after distortions and contusions, in cases of carpal-tunnel syndrome, as well as in cases of pain and overexertion; also postoperatively and in rehabilitation.

For details on materials and care see page 121.

### epX® Wrist Basic

individually packed in a transparent plastic box

Size	REF	Shipping Units
universal	22 650	1 pc. / 60 pcs.

### epX® Wrist Dynamic

individually packed in a transparent plastic box

Wrist Circumference

	Wrist Circumference	REF	Shipping Units
S	12.5 - 15.0 cm	22 660	1 pc. / 24 pcs.
M	15.0 - 16.5 cm	22 661	1 pc. / 24 pcs.
M/L	16.5 - 18.0 cm	22 662	1 pc. / 24 pcs.
L	18.0 - 20.5 cm	22 663	1 pc. / 24 pcs.
XL	20.5 - 23.0 cm	22 664	1 pc. / 24 pcs.
XXL	23.0 - 25.5 cm	22 665	1 pc. / 24 pcs.