

Cellacast® Active Bandages

Cellacast® Hook and Loop Tapes

Cast without fibreglass for focussed immobilisation. Cellacast Active allows an ideal combination of the necessary immobilisation and controlled movement. For immobilisation in primary and secondary treatment after fractures, operations, tendon injuries and overload

conditions. Cellacast Active is particularly suitable to promote mobilisation, also for children's casts.

Cellacast Active bandages can be handled fast and easily. The thin carrier material makes the bandage highly conformable to body contours. All the usual

manipulations used in casting techniques – from cast application to bivalving, windowing, and a combination with Shoe-cast – apply equally to Cellacast Active.

Cellacast Active casts are thin, light, 100% X-ray transparent, and can be easily removed and worked on with cast scissors.



Cellacast® Active

Polyester carrier material, impregnated with polyurethane resin; Bandages are rolled on a PP core and individually sealed in an laminated aluminium foil.

Cellacast® Hook and Loop Tapes

For splint fixation in focussed immobilisation, available in two versions.

Both tapes made of polyamide, dark blue, off the roll, in a dispenser.

Cellacast® Active			
individually sealed	Size	REF	Shipping Units
Creme	2.5 cm x 1.8 m	23 034	10 pcs. / 40 pcs.
	5.0 cm x 3.6 m	23 035	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	23 036	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	23 037	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	23 038	10 pcs. / 40 pcs.
Blue	5.0 cm x 3.6 m	23 040	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	23 041	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	23 042	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	23 043	10 pcs. / 40 pcs.
Orange	5.0 cm x 3.6 m	23 045	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	23 046	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	23 047	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	23 048	10 pcs. / 40 pcs.
Cellacast® Hook Tape self -adhesive			
dark blue	2.5 cm x 25 m	18 040	1 pc. / 24 pcs.
Cellacast® Loop Tape non -adhesive			
dark blue	2.5 cm x 25 m	18 041	1 pc. / 24 pcs.