

Cellacast® Soft Casting Tape

Semi-rigid casting tape for focussed immobilisation. This casting tape is suitable for making semi-rigid casts with high wearer comfort. In combination with a Cellacast Xtra splint, it allows a precise, selective increase in stability. For immobilisation in primary and secondary care of bone fractures, in postoperative care, or in the treatment of ligament injuries and conditions due to excessive physical stress. Cellacast Soft is particularly appropriate for use in the functional

mobilisation of patients and for children's casts.

Cellacast Soft tapes are quick and easy to apply. The widthwise and diagonal stretch of the material makes the application easy and convenient; its fine structure allows smooth and accurate moulding to body contours.

All common P.o.P. techniques—from application to splitting, windowing, or the combination with Cellona Shoecast—can be applied to Cellacast Soft as well.

Cellacast Soft casts are thin, lightweight, and X-ray transparent; they can be removed and trimmed to size by use of cast scissors.

Cellacast® Soft

Knitted fibreglass fabric impregnated with polyurethane resin; rolled on a reel core (PP), individually sealed in aluminium composite foil.



Cellacast® Soft Casting Tape			
	Size	REF	Shipping Units
Individually sealed, creme	2.5 cm x 1.8 m	23 721	10 pcs. / 40 pcs.
	5.0 cm x 3.6 m	23 722	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	23 723	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	23 724	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	23 725	10 pcs. / 40 pcs.
blue	5.0 cm x 3.6 m	23 726	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	23 727	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	23 728	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	23 729	10 pcs. / 40 pcs.