

Cellacast® Splints

Synthetic casting tape, padded on one side. For a clean, time saving and economical splinting technique; for functional and partially immobilising casts, e.g.

after distortions, fibular ligament rupture; in cases of tendovaginitis, epicondylitis humeri or distal fracture of the forearm.

Cellacast® Splints

7 layers of knitted fibreglass, impregnated with polyurethane resin to form the mouldable core of the splint, are integrated between fine padding felt (100% polypropylene) on the skin contact side and a film/padding cover. Cellacast splints are sealed in a moisture-proof aluminium composite foil.



Cellacast® Splints			
in a cardboard dispenser		REF	Shipping Units
folded length 3.6 m	5.0 cm width	25 039	1 box / 2 boxes
	7.5 cm width	25 040	1 box / 2 boxes
	10.0 cm width	25 041	1 box / 2 boxes
	12.5 cm width	25 042	1 box / 2 boxes
	15.0 cm width	25 043	1 box / 2 boxes
Cellacast® Pre-cut Splints			
individually sealed	7.5 x 30 cm	25 250	10 pcs. / 40 pcs.
	10.0 x 38 cm	25 251	10 pcs. / 40 pcs.
	10.0 x 76 cm	25 252	5 pcs. / 20 pcs.
	12.5 x 76 cm	25 253	5 pcs. / 20 pcs.