

# Cellacast® Xtra Casting Tape and Splints

Synthetic material for exceptionally lightweight supports with high strength; for immobilisation after fractures, operations and orthopaedic corrections; for treatment of joint and bone disorders; ideal for pro-



moting mobilisation and for bandaging child patients.

Cellacast Xtra bandages and splints are economical in material consumption. They have widthwise and diagonal stretch for easier bandaging and moulding.

All the techniques familiar from Plaster of Paris apply equally to Cellacast Xtra, e.g. for opening, splitting, windowing, combining with walking heels, hinging etc.

Cellacast Xtra bandages are lightweight, permeable to air, X-ray transparent and waterproof once hardened.



## Cellacast® Xtra

Knitted fibreglass fabric impregnated with polyurethane resin; bandages wound on a reel core (PP), splints packed in 4 layers and individually sealed in aluminium composite foil; unpadded pre-cuts for applying splints and for the reinforcement of circular casts.

Cellacast® Xtra Casting Tape			
individually sealed	Size	REF	Shipping Units
<b>Creme</b>	2.5 cm x 1.8 m	25 205	10 pcs. / 40 pcs.
	5.0 cm x 3.6 m	25 201	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	25 202	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	25 203	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	25 204	10 pcs. / 40 pcs.
<b>Blue</b>	5.0 cm x 3.6 m	25 211	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	25 212	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	25 213	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	25 214	10 pcs. / 40 pcs.
<b>Orange</b>	5.0 cm x 3.6 m	25 235	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	25 236	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	25 237	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	25 238	10 pcs. / 40 pcs.
<b>Yellow</b>	5.0 cm x 3.6 m	25 239	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	25 240	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	25 241	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	25 242	10 pcs. / 40 pcs.
<b>Green</b>	5.0 cm x 3.6 m	25 221	10 pcs. / 40 pcs.
	7.5 cm x 3.6 m	25 222	10 pcs. / 40 pcs.
	10.0 cm x 3.6 m	25 223	10 pcs. / 40 pcs.
	12.5 cm x 3.6 m	25 224	10 pcs. / 40 pcs.
Cellacast® Xtra Splints			
4-ply, individually sealed, <b>Creme</b>	5 cm x 20 cm	25 660	1. pc. / 15 pcs. / 60 pcs.
	7.5 cm x 38 cm	25 261	1. pc. / 10 pcs. / 40 pcs.
	10.0 cm x 38 cm	25 262	1. pc. / 10 pcs. / 40 pcs.
	10.0 cm x 96 cm	25 263	1. pc. / 5 pcs. / 20 pcs.
	12.5 cm x 76 cm	25 264	1. pc. / 5 pcs. / 20 pcs.