



# Cellona® Adhesive Padding Cast Edge Padding

Adhesive, very supple padding material for non-slip application under partially padded P.o.P. and synthetic casting tapes.

## Cellona® Adhesive Padding

Preferred for individually cut shapes to protect exposed support zones when applying trunk

or pelvis/leg casts. The polyacrylate adhesive is applied in strips, thus ensuring good air permeability and high skin tolerance.



## Cellona® Adhesive Cast Edge Padding

For extra padding at hard edges of casts, especially with synthetic casting tapes; preferred for forearm, upper arm and tutor casts and at all potential pressure points.

## Cellona® Adhesive Padding Cellona® Adhesive Cast Edge Padding

Nonwoven padding or air-permeable needled nonwoven: 50% polyester, 30% polypropylene, 20% viscose; polyacrylate adhesive, high skin tolerance.

<b>Cellona® Adhesive Padding</b>			
adhesive, white		REF	Shipping Units
<b>5 mm thick,</b> in a foil pouch	19 cm x 38 cm	50 810	10 pcs. / 60 pcs.
<b>2 mm thick,</b> roll form	58 cm x 1 m 28 cm x 10 m	50 811 50 813	1 box / 10 boxes 1 roll / 11 rolls
<b>Cellona® Adhesive Cast Edge Padding</b>			
white			
<b>2 mm thick,</b> roll form, in a cardboard dispenser	8 cm x 5 m	50 812	1 box / 24 boxes