



Cellona® Plaster of Paris Bandages

Ready-to-dip Plaster of Paris bandages with exceptionally fine creaminess giving enhanced mouldability; for casts applied in circular turns, for im-

mobilisation after fractures, operations and orthopaedic corrections; for treating joint and bone disorders. Cellona Plaster of Paris bandages are fully



soaked after brief immersion. The finished cast has an early weight bearing after approx. 30 minutes; it is fully set after approx. 24 hours.

Cellona® P.o.P. Bandages

Gauze bandage, porously coated with plaster compound made from extremely fine α and β semi-hydrate.

The bandages are sealed in packaging impermeable to water vapour.

Cellona® P.o.P. Bandages		REF	Shipping Units
2 m 5 bandages in a bag	6 cm	25 400	60 pcs.
	8 cm	25 401	60 pcs.
	10 cm	25 402	40 pcs.
	12 cm	25 403	40 pcs.
	15 cm	25 404	40 pcs.
	20 cm	25 405	40 pcs.
3 m	8 cm	25 411	40 pcs.
	10 cm	25 412	40 pcs.
	12 cm	25 413	40 pcs.
	15 cm	25 414	40 pcs.
	20 cm	25 415	30 pcs.
4 m	8 cm	25 421	40 pcs.
	10 cm	25 422	40 pcs.
	12 cm	25 423	30 pcs.
	15 cm	25 424	20 pcs.
	20 cm	25 425	20 pcs.
1 m, individually sealed	5 cm	24 999	10 pcs. / 200 pcs.
2.7 m individually in a bag	5 cm	91 380	12 pcs. / 96 pcs.
	7.5 cm	91 381	12 pcs. / 72 pcs.
	10 cm	91 382	12 pcs. / 48 pcs.
	12 cm	91 383	12 pcs. / 48 pcs.
	15 cm	91 384	12 pcs. / 48 pcs.
	20 cm	91 385	12 pcs. / 36 pcs.