

Cellona® Plaster of Paris Splints Wide Material

Cellona P.o.P., ready-to-dip, folded in 4 layers, or as wide material off the roll; for splints and shell-shaped plaster casts and cast sections, e.g. casts for thorax and combined thorax and

extremities, pelvis/leg casts, plaster corsets and plaster beds; also plaster casts for orthotic devices and supports. The splinting technique shortens the application time, re-

duces material consumption and is also suitable for supplementing or reinforcing casts applied in circular turns.

Cellona® P.o.P. Splints Cellona® Wide Material

Gauze bandage, porously coated with plaster compound made from extremely fine α and β semi-hydrate; 4-fold splints, zigzag folded, in PE bag impermeable to water vapour, in cardboard dispenser.



Cellona® Splints		REF	Shipping Units
4-ply			
20 m in a foil pouch and cardboard dispenser	10 cm	25 500	1 pc.
	12 cm	25 501	1 pc.
	15 cm	25 502	1 pc.
	20 cm	25 503	1 pc.
Cellona® Wide Material			
5 m roll individually in a foil pouch	40 cm width	20 301	1 pc. / 10 pcs.
	60 cm width	20 302	1 pc. / 10 pcs.
	80 cm width	20 303	1 pc. / 10 pcs.