



Cellona® Shoe

Cellona Shoe is a dressing protection shoe which can be worn over dressings either as a closed shoe or a sandal. When closed, Cellona Shoe offers reliable protection all around: the back of the foot is covered by a softly padded flap which protects foot and toes against cold and wet conditions. If it is placed onto the inside sole before putting on the shoe, Cellona

Shoe is changed into a sandal and there is free air circulation to the foot.

Cellona Shoe provides the bandaged foot with a reliable hold while walking – no matter whether the shoe is worn open or closed. The integrated hook & loop fasteners allow an individual adjustment to the foot. A buckle in the ankle area sup-



ports perfect adjustment and fit. The non-slip special profile allows safe walking with Cellona Shoe even in wet weather conditions.

Cellona Shoe can be worn both on the right and left foot. Its adjustable hook & loop fasteners make it conformable to any dressing – like wound dressings, Plaster of Paris or synthetic casts or to phlebological and lymphological dressings.

Cellona Shoe is water- and dirt-repellent and can be hand-washed at 30°C (no tumble-drying).

Cellona® Shoe			
Size		REF	Shipping Units
for children individually packed in a foil pouch	P (shoe size 25-30)	16 471	1 pc./ 12 pcs.
	XS (shoe size 31-34)	16 472	1 pc./ 12 pcs.
for adults individually packed in a foil pouch	S (shoe size 35-38)	16 473	1 pc./ 12 pcs.
	M (shoe size 39-41)	16 474	1 pc./ 12 pcs.
	L (shoe size 42-43)	16 475	1 pc./ 12 pcs.
	XL (shoe size 44-47)	16 476	1 pc./ 12 pcs.