

Cellona® Shoecast • Walking Heel

For physiologically correct movement with a walking cast and for protection against moisture and dirt, L&R offer two walking soles:



Cellona® Shoecast

Walking sole with a high degree of lengthwise and widthwise stability plus maximum ground contact. The anatomical shape permits natural heel-to-toe movement. The highly abrasion-resistant sole profile prevents slipping, while the foamed material cushions every step. A built-up heel cup and cross-over straps with infinitely adjustable hook & loop fastener ensure a secure fit.

Walking sole: 100% foamed polyurethane; straps: 100% polyamide.



Cellona® Walking Heel

Profiled walking heel cast from light grey synthetic rubber; lightweight, anti-slip and shock-absorbing; grooves and lugs for attachment either immediately or after application of P.o.P. or synthetic lower limb casts.



Cellona® Shoecast			
		REF	Shipping Units
individually in a foil pouch			
for children	Size 0, left	50 858	1 pc. / 8 pcs.
	Size 0, right	50 859	1 pc. / 8 pcs.
for adults	Size 1, small left	50 860	1 pc. / 10 pcs.
	Size 1, small right	50 861	1 pc. / 10 pcs.
	Size 2, medium left	50 862	1 pc. / 10 pcs.
	Size 2, medium right	50 863	1 pc. / 10 pcs.
	Size 3, large left	50 864	1 pc. / 8 pcs.
	Size 3, large right	50 865	1 pc. / 8 pcs.
	Size 4, extra large left	50 869	1 pc. / 6 pcs.
	Size 4, extra large right	50 870	1 pc. / 6 pcs.
Replacement Strap with Hook & Loop Fastener	40 cm (for size 0)	50 868	10 pcs. / 80 pcs.
	90 cm (for sizes 1 and 2)	50 873	10 pcs. / 100 pcs.
	100 cm (for sizes 3 and 4)	50 874	10 pcs. / 100 pcs.
Cellona® Walking Heel			
small		50 820	1 pc. / 20 pcs.
medium		50 821	1 pc. / 20 pcs.
large		50 822	1 pc. / 10 pcs.