



Cellona® Synthetic Undercast Padding

Undercast padding that retains its softness and elasticity to protect exposed bones and nerves under P.o.P. and synthetic cast-

ing tapes. The fine fibrous web is strong enough to permit application in circular turns, but is easy to tear and provides a

smooth surface. The material's self-adhesion facilitates application by simply smoothing down.

Cellona synthetic undercast padding is air-permeable and compensates for changes in temperature. The fibres do not absorb any moisture, the casting tape hardens faster, and the good padding effect is retained with long-term wear.



Cellona® Synthetic Undercast Padding

100% polyester fibres, no optical brighteners, in various sizes in roll form, with paper wrapper; sterilisable (A 121°C; or other methods).

Cellona synthetic undercast padding is also used for padding splints, zinc paste bandages (Unna's Boot) and adhesive bandages. For wound treatment – in combination with textile dressings or ointment tulle - Cellona padding acts as an air-permeable cushion.

Cellona® Synthetic Undercast Padding

	REF	Shipping Units
6 cm x 3 m	10 686	6 rolls / 60 rolls
10 cm x 3 m	10 687	4 rolls / 40 rolls
40 cm x 10 m	10 697	1 roll / 4 rolls
4 cm x 3 m	10 690	108 rolls
6 cm x 3 m	10 692	72 rolls
10 cm x 3 m	10 694	48 rolls
15 cm x 3 m	10 695	36 rolls
20 cm x 3 m	10 696	24 rolls
5 cm x 2.7 m	92 840	96 rolls
7.5 cm x 2.7 m	92 841	72 rolls
10 cm x 2.7 m	92 842	48 rolls
15 cm x 2.7 m	92 843	36 rolls
20 cm x 2.7 m	92 844	24 rolls