



Cellona® Xtra P.o.P. Bandages and Splints

The main difference between Cellona and Cellona Xtra Plaster of Paris is the shorter setting time of Cellona Xtra.

Simple dip technique: with easy-soak through properties and low plaster loss.

Enhanced mouldability: extra plaster with extra creaminess. Guaranteed 2.5 minutes setting time.

Windowing and bivalving after 10 minutes demonstrate excellent early strength.

Partial weight bearing after half an hour, with assistance of walking aid.

Full weight bearing after 24 hours.

Cellona Xtra Plaster of Paris bandages and splints are indicated for the immobilisation of fractures, postoperatively after orthopaedic surgical procedures, and in the treatment of bone and joint diseases.



Cellona® Xtra Plaster of Paris Bandages

	Size	REF	Shipping Units
individually sealed	5 cm x 2.7 m	91 180	1 box / 60 pcs.
	7.5 cm x 2.7 m	91 181	1 box / 60 pcs.
	10 cm x 2.7 m	91 182	1 box / 60 pcs.
	15 cm x 2.7 m	91 184	1 box / 36 pcs.
	20 cm x 2.7 m	91 185	1 box / 36 pcs.
	6 cm x 3 m	25 050	1 box / 60 pcs.
	8 cm x 3 m	25 051	1 box / 60 pcs.
	10 cm x 3 m	25 052	1 box / 60 pcs.
	12 cm x 3 m	25 053	1 box / 48 pcs.
	15 cm x 3 m	25 054	1 box / 36 pcs.
	20 cm x 3 m	25 055	1 box / 30 pcs.

Cellona® Xtra Plaster of Paris Splints, 4-ply

each splint sealed in a protective PE bag and packed to cardboard dispenser	10 cm x 20 m	25 056	1 box / 1 splint
	12 cm x 20 m	25 057	1 box / 1 splint
	15 cm x 20 m	25 058	1 box / 1 splint
	20 cm x 20 m	25 059	1 box / 1 splint