

Cotton Padding Padding Bandages

Padding in various combinations of fibres for partial, all-over or circular application under P.o.P. and synthetic casts.

Cotton Padding

Approx. 80% cotton, not degreased, approx. 20% viscose, on a roll, with paper layer.

Padding Bandages

Padding for medical purposes DAB (German Pharmacopoeia), 50% bleached cotton, 50% viscose, in rolls with paper wrapper.



Cotton Padding			
off-white	Size	REF	Shipping Units
individually in a pouch	40 cm x 2.40 m	10 603	1 pack / 20 packs
Padding Bandages			
Padding for medical purposes DAB	10 cm x 5 m 15 cm x 5 m	10 640 10 641	1 box (75 rolls) 1 box (50 rolls)