

Haftan® Crepe Paper Bandages

Haftan®

Extra-thin, supple bandage made from white polyurethane foam which does not slip and creates a smooth, even surface, holding padding material in place under casting tapes, es-

pecially in combination with synthetic material.

The fine, open-pored surface of Haftan allows the individual turns of the bandage to adhere well to each other, thus preventing slippage even on the skin.

The PU foam reduces any penetration of moisture.

Haftan is also used for joint supports (sports), combined with adhesive bandages or tape.



Crepe Paper Bandage

Fine, crimped paper bandage, lengthwise slightly stretchable. Allows crease-free application; for protecting the skin under unpadded plaster casts and for holding padding material in place. The bandage retains its strength even when wet.

100% cellulose

Haftan®		REF	Shipping Units
white, thickness: 0.8 mm	7 cm x 27.5 m	22 315	6 pcs. / 60 pcs.
	10 cm x 27.5 m	22 316	6 pcs. / 60 pcs.
Crepe Paper Bandages			
4 m	4 cm	14 250	20 pcs. / 200 pcs.
	6 cm	14 251	20 pcs. / 200 pcs.
	8 cm	14 252	20 pcs. / 200 pcs.
	10 cm	14 253	20 pcs. / 200 pcs.
	12 cm	14 254	20 pcs. / 200 pcs.