

tg[®] Tubular Bandage

A fine, seamlessly knitted tube with widthwise extensibility for use as crease-free undercast protection – from the little finger through to body-size P.o.P. or synthetic casts; also suitable to hold padding material in place, and for dermatological indications.

tg[®] Tubular Bandage

67% bleached cotton, 33% viscose, sterilisable (A 121°C).



tg tubular bandage is also used for dermatological indications, over ointment dressings or zinc paste bandages or as undercast protection for P.o.P. or synthetic casts.

tg[®] Tubular Bandage				REF	Shipping Units
white, individually in a cardboard dispenser					
5 m roll	size 1	(1.4 cm)	for finger and toe	24 020	1 roll / 80 rolls
	size 2	(2.3 cm)	for fingers and toes	24 021	1 roll / 50 rolls
	size 3	(3 cm)	for several fingers, child's arm	24 022	1 roll / 40 rolls
	size 5	(5.5 cm)	for arm, lower leg, child's leg	24 023	1 roll / 40 rolls
	size 6	(6.5 cm)	for arm, leg (medium)	24 024	1 roll / 40 rolls
	size 7	(7 cm)	for arm (large), leg, child's head	24 025	1 roll / 40 rolls
	20 m roll	size 1	(1.4 cm)	for finger and toe	24 000
size 2		(2.3 cm)	for fingers and toes	24 001	1 roll / 20 rolls
size 3		(3 cm)	for several fingers, child's arm	24 002	1 roll / 10 rolls
size 5		(5.5 cm)	for arm, lower leg, child's leg	24 003	1 roll / 10 rolls
size 6		(6.5 cm)	for arm, leg (medium)	24 004	1 roll / 10 rolls
size 7		(7 cm)	for arm (large), leg, child's head	24 005	1 roll / 18 rolls
size 9		(8.5 cm)	for thigh (large), head, armpit	24 006	1 roll / 12 rolls
size 12		(12 cm)	for thigh (extra large), head (large), child's trunk	24 007	1 roll / 12 rolls
10 m roll	size K1	(16 cm)	for trunk (small) < M	24 030	1 roll / 10 rolls
	size K2	(21 cm)	for trunk (large) > L	24 031	1 roll / 10 rolls