

tg® Vest

Short-sleeved vest, soft, supple and skin-friendly, as an easy-to-apply and time-saving under-cast protection for P.o.P. and synthetic casting tapes. The elastic knitted material has widthwise extensibility, allow-

ing for intermediate sizes. If two vests are applied – and a special casting technique is used – the vest that is in contact with the skin can be changed, an advantage for body hygiene under long-term casts.

tg vests are also used in dermatology to cover large areas of affected skin.

100% cotton, not degraded or optically brightened, sterilisable (A 121°C).



tg® Vest

	REF	Shipping Units
individually in a cardboard box		
small (ladies' size 36 - 44, men 42 - 50)	24 044	1 pc. / 14 pcs.
large (ladies' up to size 46, men up to size 56)	24 045	1 pc. / 9 pcs.