


Raucodrape®

High Performance for your OR

The new Raucodrape OR drapes meet the highest standards in terms of EN 13795.

The safety features of Raucodrape®

- safe barrier function against liquid penetration
- safe barrier function against bacterial penetration
- superior tear resistance - even when wet
- high cleanliness both microbial and particulate
- extremely low linting




Raucodrape® with even higher tear-resistance and good absorption capacity

* PE: polyethylene, PP: polypropylene

Tear resistance


Raucodrape has excellent tear resistance. The improvement achieved in the wet material is outstanding.



all values given in Newton N; tested in accordance with ISO 9073-3
*as compared to the similar previous version of the material

Water absorption capacity

Raucodrape can absorb 250 ml of water per m².



tested in accordance with ISO 9073-6

Liquid penetration

Raucodrape is liquid-proof. Even under a water column of 125 cm, it still forms an impenetrable barrier.



tested in accordance with EN 20811