


Raucodrape® OR Drape System

Cardiology/Angiography

Raucodrape® Cardiovascular Set, sterile			
Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 380	7 sets
4 towels	30 x 40 cm		
2 instrument table cover	80 x 145 cm		
2 adhesive tapes	9 x 50 cm		
2 foot protectors	37.5 x 24 cm		
1 OR drape sheet	150 x 180 cm		
1 split sheet bilateral	<u>175</u> x 260 cm		
2-layered Plus, split, with integrated tube retainer	20 x 100 cm		
1 cardiovascular drape	260/200 x 335 cm		
2-/3-layered, aperture, with integrated incision film, with 3-chamber instrument pouches on both sides, with integrated tube retainer and a cord organiser	33 x 38 cm		
Raucodrape® Angiography Drape I, sterile			
1 wrapping paper	80 x 100 cm	33 381	11 sets / 22 sets
4 towel	30 x 40 cm		
1 adhesive tape	9 x 25 cm		
1 angiography drape, 3-layered	150 x 240 cm		
2 apertures (circular, off-centre)	12 cm		
Raucodrape® Angiography Drape II, sterile			
1 wrapping paper	80 x 100 cm	33 382	8 sets / 16 sets
1 towel	30 x 40 cm		
1 adhesive tape	9 x 25 cm		
1 angiography drape, 3-layered, PE film on both sides of drape,	227 x 330 cm		
2 apertures (off-centre) with fenestrated incision film	12 x 12 cm 8 cm		
 adhesive side			