

Raucodrape® OR Drape System

Neurosurgery/ENT/CMF/ Ophthalmology

Raucodrape® Kranio-Set, sterile			
Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 370	14 sets
2 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
3 adhesive drapes (PE film)	50 x <u>50</u> cm		
1 craniotomy drape, aperture,	225 x 280 cm		
with integrated incision film, integrated fluid collection pouch with lateral gutters, sieve, and exit port, with integrated tube retainers	19 x 25 cm		
Raucodrape® LAP/Spinal Drape, sterile			
LAP/Spinal drape, 3-layered aperture (oval, lengthwise, off-centre, self-adhesive)	225 x 320 cm 11 x 22 cm	33 051	8 pcs. / 16 pcs.
Raucodrape® Head Set I, sterile			
1 table cover	150 x 190 cm	33 371	16 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 OR drape sheet, 2-layered	75 x 90 cm		
1 split sheet, 2-layered, split	150 x 250 cm 10 x 80 cm		
1 adhesive drape, 2-layered	75 x <u>90</u> cm		
Raucodrape® Head-Set II, sterile			
1 table cover	150 x 190 cm	33 372	13 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 adhesive tape	9 x 50 cm		
1 turban drape (2-piece) fluid-repellent nonwoven material,	<u>130</u> x 100 cm 225 x 300 cm		
2-layered,	14 x 80 cm		
1 split sheet, 2-layered Plus, split (rectangular), with integrated tube retainers			
— adhesive side			



Raucodrape® OR Drape System

Neurosurgery/ENT/CMF/ Ophthalmology

Raucodrape® Variable Split Sheet Set, sterile			
Contents	Size	REF	Shipping Units
1 table cover	150 x 190 cm	33 313	15 sets
4 towels	30 x 40 cm		
1 instrument table cover	80 x 145 cm		
1 split sheet, 2-layered Plus, split	225 x 240 cm 20 x 20 cm		
1 split sheet, 2-layered Plus, split	150 x 260 cm 20 x 60 cm		
Raucodrape® Split Sheet, sterile			
Split Sheet, 2-layered split	150 x 250 cm 10 x 80 cm	33 073	14 pcs. / 28 pcs.
Raucodrape® Fenestrated Face Drape, sterile			
Fenestrated face drape, 2-layered, aperture, with integrated incision film, (oval, widthwise, off-centre)	150 x 240 cm 8 x 19 cm	33 054	15 pcs. / 30 pcs.
Raucodrape® Fenestrated Ear Drape, sterile			
Fenestrated ear drape, 2-layered, aperture (oval, lengthwise, off-centre, self-adhesive)	175 x 240 cm 5 x 7.5 cm	33 050	11 pcs. / 22 pcs.
Raucodrape® Eye Drape, sterile			
Eye drape (fluid-repellent nonwoven material), central aperture, with integrated incision film, integrated fluid collection pouch	150 x 150 cm 7 x 10 cm	33 056	17 pcs. / 34 pcs.
— adhesive side			