



# Sentinex<sup>®</sup> Surgical Face Masks

## Sentinex<sup>®</sup> Surgical Caps



Sentinex lite



Sentinex CLASSIC



Sentinex SuperSoft



Sentinex Extra Touch



Sentinex Safety Shield with anti-fog visor



Sentinex Ultra



Sentinex AntiFog



Sentinex Soft



Sentinex Fancy



Sentinex Astro



Sentinex Astro Extra



Sentinex Color



Sentinelx Nova



Sentinelx Trend



Sentinelx Flex



Sentinelx Easy

Sentinelx Surgical Face Masks in different versions, made as if “to measure”, e.g. with optimum wearer comfort and proven BFE value\* for efficient prophylaxis against infections for patients and staff, in the OR and doctor’s surgery, in intensive care units and in areas with high infection tendency. The extremely long nose bow is made of plastic covered steel

and is integrated into soft, skin-friendly seams and offers good fit around the nose. Filtration material, inner and outer layers and the ties are ultrasound sealed – a hygienically perfect and high quality finish. All Sentinelx Surgical Face Masks and Sentinelx Surgical Caps are free of latex and fiberglass.

### Sentinelx® Surgical Face Masks

	REF	Shipping Units
Sentinelx Lite, blue	11 983	50 pcs. / 300 pcs.
Sentinelx Lite, green	11 991	50 pcs. / 300 pcs.
Sentinelx Lite, green, with ear loops, latex-free	11 992	50 pcs. / 500 pcs.
Sentinelx CLASSIC, green	11 996	50 pcs. / 300 pcs.
Sentinelx SuperSoft	11 997	50 pcs. / 300 pcs.
Sentinelx Extra Touch	11 998	50 pcs. / 300 pcs.
Sentinelx Safety Shield	11 989	25 pcs. / 100 pcs.
Sentinelx AntiFog	11 986	50 pcs. / 300 pcs.
Sentinelx Ultra	11 999	50 pcs. / 300 pcs.

### Sentinelx® Surgical Caps

Sentinelx Soft	Beret cap, green, breathable	17 420	100 pcs. / 600 pcs.
Sentinelx Astro	Astro cap, green	17 421	100 pcs. / 600 pcs.
Sentinelx Astro Extra	Astro cap with sweatband, green	17 422	100 pcs. / 600 pcs.
Sentinelx Fancy	Nurse’s cap with center seam, green	17 423	100 pcs. / 600 pcs.
Sentinelx Color	Nurse’s cap with center seam, colored	17 430	100 pcs. / 600 pcs.
Sentinelx Nova	Nurse’s cap with strings, green	17 424	100 pcs. / 600 pcs.
Sentinelx Trend	Nurse’s cap, long, green	17 425	100 pcs. / 600 pcs.
Sentinelx Flex	Surgeon’ cap, green	17 426	100 pcs. / 600 pcs.
Sentinelx Easy, green	Beret cap, green	17 427	100 pcs. / 600 pcs.
Sentinelx Easy, blue	Beret cap, blue	17 428	100 pcs. / 600 pcs.
Sentinelx Easy, white	Beret cap, white	17 429	100 pcs. / 600 pcs.