



Sentinex® Surgical Gowns • OR Stockings

Sentinex Surgical gowns by L&R represent a clearly structured concept:

- Spunlace and SMMS
- Basic, Standard, Special versions
- With and without sterile wrapping ("direct")

Sentinex OR gowns offer high wearer comfort and protection for users and patients.




Two materials: Spunlace and SMMS

With the introduction of EN 13 795 uniform requirements have been defined for OR gowns. The central requirement of this standard is the compliance with a barrier function against fluids and germs in order to ensure protection against infections, for the patient as well as for the user.


Both materials fulfill the requirements of the standard and are suitable for all types of surgical procedures – if reinforced through additional reinforcements in the critical areas.

Spunlace

- excellent wearing comfort
- skin-friendly and breathable
- ideally suited for physical demanding operations
- best possible temperature exchange




Surgical gown with sterile wrapping



- ideally suited for sterile preparation
- with sterile wrapping
- 2 towels

SMMS

- light
- protects reliably against bacteria and fluids
- guarantees the circulation of air, heat, and water vapour
- reduces thermo-physiological stress



Surgical gown "direct"



- for immediate sterile handover
- for fast removal
- without sterile wrapping and towels

Sentinex[®] Surgical Gown Basic

Gown with tie-closure for small operations and interventions, with 3 ties for a safe closure in the back.



Sentinex[®] Surgical Gown Standard

Wraparound surgical gown with a sterile back and a specially wide overlapping area. Ideally suited for operations with small amounts of fluids.



Sentinex[®] Surgical Gown Special

Wraparound surgical gown on the basis of the Standard version, with additional impermeable reinforcement in the front and sleeve area. Ideally suited for operations with a large amount of fluids.



Sentinex[®] Surgical Gown Spunlace, sterile

Contents	Size	REF	Shipping Units
Standard with sterile wrapping	120 cm	19 300	44 pcs.
	130 cm	19 301	40 pcs.
	150 cm	19 302	36 pcs.
	170 cm	19 303	32 pcs.
	150 cm large	19 304	32 pcs.
Special with sterile wrapping	120 cm	19 310	40 pcs.
	130 cm	19 311	36 pcs.
	150 cm	19 312	32 pcs.
	170 cm	19 313	28 pcs.
	150 cm large	19 314	28 pcs.
Standard "direct"	130 cm	19 321	40 pcs.
	150 cm	19 322	36 pcs.
Special "direct"	130 cm	19 331	36 pcs.
	150 cm	19 332	32 pcs.

Sentinex[®] Surgical Gowns SMMS, sterile

Standard with sterile wrapping	120 cm	19 340	52 pcs.
	130 cm	19 341	48 pcs.
	150 cm	19 342	44 pcs.
	170 cm	19 343	40 pcs.
	150 cm large	19 344	40 pcs.
Special with sterile wrapping	120 cm	19 350	48 pcs.
	130 cm	19 351	44 pcs.
	150 cm	19 352	40 pcs.
	170 cm	19 353	36 pcs.
	150 cm large	19 354	36 pcs.
Standard "direct"	120 cm	19 360	52 pcs.
	130 cm	19 361	48 pcs.
	150 cm	19 362	44 pcs.
Special "direct"	120 cm	19 380	48 pcs.
	130 cm	19 381	44 pcs.
	150 cm	19 382	40 pcs.

Sentinex[®] Surgical Stockings

Sentinex Surgical Stocking	20 225	200 pcs. / 200 pcs.
Sentinex Surgical Sock	20 224	200 pcs. / 400 pcs.