

# Setpack® OR Dressing Materials

## Flexible combination of standardised components

Thorough preparation takes its time – also in emergency situations. What is required is immediate accessibility and safety.

The Setpack Standard System offers ready-for-use sets which contain the **surgical dressing materials** most commonly used in operations – sterile and ready at any time.

Ready-for-use Standard Sets contain **wound dressings, dissecting swabs, gauze balls or laparotomy swabs** in quantities as they are needed. All Setpack products are provided with an X-ray detectable marker (min. 60% barium sulfate). Should a case of doubt arise, an X-ray control will give immediate



clarification. Depending on the requirements for your specific operations, you can make your individual choice. Determined quantities offer maximum economy.



## Patented Double Counting Card System

- dual check system by having each counting card both inside the pack and in the peel bag
- risk-free product handover during the surgery
- clear documentation by the use of self-adhesive counting cards

## Special covering of the inside pack

- a handy flap makes it easy to open the package
- fast product removal through special packaging design

## Patented packaging box for gauze balls

- increased safety by optimized counting control before and after the surgery
- special developed packaging design for handsome removal of dissecting swabs



# Setpack® OR Dressing Materials Standard Sets

<b>Setpack® Gauze Dressings, X-ray detectable, sterile</b>						
		Size	REF	Shipping Units		
Ready-for-use dressings made of absorbent gauze type 17 (17 threads/cm <sup>2</sup> ) for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card.	<b>Sets of 10 pcs.</b>	12-ply	7.5 x 7.5 cm	15 114	28 / 84 sets	
			10 x 10 cm	15 100	20 / 60 sets	
			10 x 12.5 cm	15 116	20 / 60 sets	
			10 x 20 cm	15 102	16 / 48 sets	
	16-ply	5 x 5 cm	15 110	34 / 102 sets		
		7.5 x 7.5 cm	15 106	24 / 72 sets		
		10 x 10 cm	15 101	20 / 60 sets		
		10 x 12.5 cm	15 111	18 / 54 sets		
Absorbent gauze type 20 (20 threads/cm <sup>2</sup> ) for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card.	<b>Sets of 10 pcs.</b>	32-ply	10 x 10 cm	15 104	12 / 36 sets	
Type 17 (17 threads/cm <sup>2</sup> ) for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card.	<b>Sets of 20 pcs.</b>	12-ply	7.5 x 7.5 cm	15 119	20 / 60 sets	
			10 x 10 cm	15 103	18 / 54 sets	
			10 x 12.5 cm	15 109	10 / 30 sets	
			10 x 20 cm	15 108	8 / 24 sets	
	16-ply	5 x 5 cm	15 117	26 / 78 sets		
		7.5 x 7.5 cm	15 118	16 / 48 sets		
		10 x 10 cm	15 105	14 / 42 sets		
		10 x 12.5 cm	15 112	8 / 24 sets		
<b>Setpack® Gauze Balls, X-ray detectable, sterile</b>						
Ready-for-use gauze balls made of absorbent gauze type 20 (20 threads/cm <sup>2</sup> ) with elastic, for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card	<b>Sets of 5 pcs.</b>	fist size	7	22 800	16 / 48 sets	
	<b>Sets of 10 pcs.</b>	walnut size	2	22 801	72 / 144 sets	
			plum size	3	22 802	38 / 114 sets
			egg size	4	22 803	26 / 78 sets
			extra large	6	22 804	18 / 54 sets
	<b>Sets of 15 pcs.</b>	plum size	3	22 805	30 / 90 sets	
			egg size	4	22 806	20 / 60 sets
			extra large	6	22 807	12 / 36 sets
	<b>Sets of 20 pcs.</b>	plum size	3	22 808	26 / 78 sets	
			egg size	4	22 809	16 / 48 sets
			extra large	6	22 810	6 / 18 sets
<b>Setpack® Dissecting Swabs, sterile</b>						
Ready-for-use dissecting swabs made of absorbent gauze type 24 (24 threads/cm <sup>2</sup> ) for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card.	<b>Sets of 10 pcs.</b>	small	2	12 780	40 / 120 sets	
		medium	3	12 781	40 / 120 sets	
		large	4	12 782	40 / 120 sets	
	<b>Sets of 20 pcs.</b>	small	2	12 783	40 / 120 sets	
		medium	3	12 784	40 / 120 sets	
		large	4	12 785	40 / 120 sets	

# Setpack® OR Dressing Materials

## Setpack® Toptex® lite Laparotomy Swabs, X-ray detectable

Ready-for-use laparotomy swabs made of absorbent gauze type 20 (20 threads/cm<sup>2</sup>) EN 14079, prewashed, for use in the OR, with X-ray detectable chip and self-adhesive, double, separate counting card.



### Toptex® lite Laparotomy Swabs, X-ray detectable, sterile

Ready-for-use laparotomy swabs made of absorbent gauze type 20 (20 threads/cm<sup>2</sup>), prewashed, for use in the OR, with X-ray detectable chip and double, separated, adhesive counter card.

			Size	REF	Shipping Units	
White	set of 2 pcs.	4 layers	40 x 40 cm	15 004	14 / 42 sets	
	set of 5 pcs.	4 layers		15 001	6 / 18 sets	
	set of 2 pcs.	6 layers		15 012	12 / 36 sets	
	set of 5 pcs.	6 layers	20 x 30 cm	15 003	6 / 18 sets	
	set of 2 pcs.	4 layers		15 011	28 / 84 sets	
	set of 5 pcs.	4 layers		15 007	14 / 42 sets	
	set of 2 pcs.	4 layers		50 x 60 cm	15 014	8 / 24 sets
	set of 2 pcs.	4 layers			15 010	26 / 78 sets
	set of 2 pcs.	4 layers			40 x 40 cm	15 006
	Green	set of 5 pcs.	4 layers	15 002		6 / 18 sets
set of 2 pcs.		6 layers	15 008	12 / 36 sets		
set of 5 pcs.		6 layers	15 000	6 / 18 sets		
set of 2 pcs.		4 layers	20 x 30 cm	15 017		28 / 84 sets
set of 5 pcs.		4 layers		15 005		14 / 42 sets
set of 2 pcs.		4 layers		15 013		8 / 24 sets
set of 2 pcs.		6 layers	50 x 60 cm	15 018		6 / 18 sets
set of 2 pcs.		4 layers		15 009		26 / 78 sets
Blue		set of 2 pcs.		4 layers	40 x 40 cm	15 016
	set of 5 pcs.	4 layers	15 015	6 / 18 sets		

### Toptex® lite Disposable Laparotomy Swabs, X-ray detectable, non-sterile

made of absorbent gauze Ph. Eur. type 20 (20 threads/cm<sup>2</sup>) with 30 cm loop, prewashed, with X-ray detectable chip.

4 layers	White	8 x 90 cm	10 921	100 / 600 pcs.	
		20 x 30 cm	10 922		
		40 x 40 cm	10 923		
		50 x 60 cm	10 925		
	Green	8 x 90 cm	10 916	100 / 600 pcs.	
		20 x 30 cm	10 917		
		40 x 40 cm	10 918		
		50 x 60 cm	10 920		
	Blue	40 x 40 cm	10 915	90 / 180 pcs.	
		White	40 x 40 cm	10 924	75 / 150 pcs.
	6 layers		Green	40 x 40 cm	10 919