



# Setpack® OR Dressing Materials Standard Sets

<b>Setpack® Gauze Dressings, X-ray detectable, sterile</b>					
		Size	REF	Shipping Units	
Ready-for-use dressings made of absorbent gauze type 17 (17 threads/cm <sup>2</sup> ) for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card.	<b>Sets of 10 pcs.</b>	12-ply	7.5 x 7.5 cm	15 114	28 / 84 sets
			10 x 10 cm	15 100	20 / 60 sets
			10 x 12.5 cm	15 116	20 / 60 sets
			10 x 20 cm	15 102	16 / 48 sets
	16-ply	5 x 5 cm	15 110	34 / 102 sets	
		7.5 x 7.5 cm	15 106	24 / 72 sets	
		10 x 10 cm	15 101	20 / 60 sets	
		10 x 12.5 cm	15 111	18 / 54 sets	
Absorbent gauze type 20 (20 threads/cm <sup>2</sup> ) for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card.	<b>Sets of 10 pcs.</b>	32-ply	10 x 10 cm	15 104	12 / 36 sets
Type 17 (17 threads/cm <sup>2</sup> ) for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card.	<b>Sets of 20 pcs.</b>	12-ply	7.5 x 7.5 cm	15 119	20 / 60 sets
			10 x 10 cm	15 103	18 / 54 sets
			10 x 12.5 cm	15 109	10 / 30 sets
			10 x 20 cm	15 108	8 / 24 sets
	16-ply	5 x 5 cm	15 117	26 / 78 sets	
		7.5 x 7.5 cm	15 118	16 / 48 sets	
		10 x 10 cm	15 105	14 / 42 sets	
		10 x 12.5 cm	15 112	8 / 24 sets	
<b>Setpack® Gauze Balls, X-ray detectable, sterile</b>					
Ready-for-use gauze balls made of absorbent gauze type 20 (20 threads/cm <sup>2</sup> ) with elastic, for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card	<b>Sets of 5 pcs.</b>	fist size	7	22 800	16 / 48 sets
	<b>Sets of 10 pcs.</b>	walnut size	2	22 801	72 / 144 sets
			3	22 802	38 / 114 sets
			4	22 803	26 / 78 sets
			6	22 804	18 / 54 sets
	<b>Sets of 15 pcs.</b>	plum size	3	22 805	30 / 90 sets
			4	22 806	20 / 60 sets
			6	22 807	12 / 36 sets
	<b>Sets of 20 pcs.</b>	plum size	3	22 808	26 / 78 sets
			4	22 809	16 / 48 sets
			6	22 810	6 / 18 sets
<b>Setpack® Dissecting Swabs, sterile</b>					
Ready-for-use dissecting swabs made of absorbent gauze type 24 (24 threads/cm <sup>2</sup> ) for use in the OR, with X-ray detectable thread and double, separated, adhesive counter card.	<b>Sets of 10 pcs.</b>	small	2	12 780	40 / 120 sets
		medium	3	12 781	40 / 120 sets
		large	4	12 782	40 / 120 sets
	<b>Sets of 20 pcs.</b>	small	2	12 783	40 / 120 sets
		medium	3	12 784	40 / 120 sets
		large	4	12 785	40 / 120 sets