



Curatest® Curatest® F

Curatest® Nonwoven Patch Test Strips

The patch test strip with 10 test discs is as simple to apply as a first-aid dressing. The detachable side edges can be written on to assist diagnosis.

Nonwoven carrier material: 100% polypropylene
Frame of test discs: polyester film
Absorbent test discs: 60% cotton, 40% viscose
Polyacrylate adhesive: free of colophonium and colophonium derivatives, hypoallergenic.

Curatest® F Film Patch Test Strips

This patch test strip with 5 or 10 test discs has a semi-permeable, transparent carrier film that is permeable to water vapour but waterproof. It adheres reliably, like a second skin, even when taking a shower. Easy application of the test substances to the discs and of the test strip to the skin.

Carrier film: 100% polyurethane.
Frame of test discs: polyester film.
Absorbent test discs: 60% cotton, 40% viscose.
Polyacrylate adhesive: free of colophonium and colophonium derivatives, hypoallergenic.



Curatest® Nonwoven Patch Test Strips

Size	REF	Shipping Units
each strip has 10 test discs 7.5 x 12.5 cm	30 054	50 pcs. / 1200 pcs.

Curatest® F Film Patch Test Strips

each strip has 10 test discs 7.5 x 12.5 cm	30 062	50 pcs. / 1200 pcs.
each strip has 5 test discs 3.2 x 12.5 cm	30 061	100 pcs. / 2000 pcs.