

Impregnated Ribbon Gauze

Paraffin Ribbon Gauze

Vaseline-impregnated ribbon gauze with selvages for packing cavity wounds, after surgical interventions and for application in the fields of ear, nose and throat treatment and orthodontic surgery.

The impregnation with white vaseline is dosed in such a way that no vaseline residue remains in the wound after removal.



Carrier material: Ribbon gauze type 24 (24 threads/cm²) with selvages, 100% cotton

Ointment: White vaseline (acc. to Ph. Eur.)

Clauden® Nasal Ribbon Gauze

The impregnation with Clauden®* powder and other substances prevents adhesion of the ribbon gauze to the wound tissue.

Clauden® Nasal Ribbon Gauze consists of:

- Clauden powder
- Clioquinol
- Liquid paraffin
- White vaseline
- Lanolin
- Talcum
- Ribbon gauze (cotton)



100% cotton (carrier material) ointment basis in laminated aluminium pouch.

* Trademark, under licence from DAIICHI SANKYO EUROPE GmbH

Opraclean® Ribbon Gauze

Iodoform impregnated ribbon gauze, antiseptic, can be applied to the wound precisely in the required dosis.

Practical dispensers make these dressings easy to handle in everyday situations in hospitals and doctor's surgeries; at the same time they satisfy hygiene requirements.



Paraffin Ribbon Gauze

	REF	Shipping Units
2 cm x 1 m, individually packed in aluminium pouch	13 514	10 pcs. / 300 pcs.
2 cm x 10 m, individually packed in plastic jar	13 515	1 pc. / 96 pcs.

Clauden® Nasal Ribbon Gauze

	REF	Shipping Units
2 cm x 1 m, in laminated aluminium pouch	21 000	1 pc. / 240 pcs.

Opraclean® Ribbon Gauze

	REF	Shipping Units
1 cm x 5 m, rolled, in dispenser	13 580	1 pc. / 216 pcs.
2 cm x 5 m, rolled, in dispenser	13 581	1 pc. / 120 pcs.
5 cm x 5 m, rolled, in dispenser	13 582	1 pc. / 60 pcs.