



Pro-ophta® Eye Cushion Pad

Pro-ophta® Eye Dressing Pad

Pro-ophta® Eye Cushion Pad

Eye cushion with sealed edges, soft and lint-free for use in the treatment of eye injuries, post-operative eye care or after application of ointments.

Interior: absorbent, soft wadding, 60% polypropylene, 40% cotton; exterior: soft nonwoven casing made of 100% polypropylene.

Pro-ophta® Eye Dressing Pad

Oval eye dressing pad, particularly soft with surgical gauze on both sides. For treatment of external eye injuries, postoperative eye care or after the application of ointments.

Interior: eye-grade cotton wool DIN; exterior: surgical gauze casing, 28 threads, 100% cotton; oval shape, non sterile with paper liner.



Pro-ophta® Eye Cushion Pad			
		REF	Shipping Units
sterile , individually sealed	6.2 x 7.2 cm	13 043	25 pcs. / 250 pcs.
non-sterile in a cardboard box	6.2 x 7.2 cm	13 041	50 pcs. / 500 pcs.
Pro-ophta® Eye Dressing Pad			
sterile individually sealed	5.5 x 7.5 cm	13 019	5 pcs. / 200 pcs.
non-sterile in a cardboard box	5.5 x 7.5 cm	13 020 13 021	5 pcs. / 250 pcs. 50 pcs. / 500 pcs.