



Pro-ophta® Eye Dressings

Pro-ophta® Junior

The skin suffers when adhesive dressings have to be changed frequently. Pro-ophta Junior helps protect children's delicate skin by means of an extremely thin nonwoven and a porously

applied adhesive. The dressing allows plenty of air to reach the skin and is easy to remove without leaving any residues, without pulling or tearing. Even children with sensitive skin therefore experience no discomfort with Pro-ophta Junior.

The occlusive cotton centre section keeps the eye in complete darkness. With its shape, specially designed for a child's eye, Pro-ophta Junior offers a perfect fit. Whether it is worn on the left or right eye, with or without spectacles, Pro-ophta Junior adapts optimally to the shape of the eye, and the wearer scarcely notices it thanks to the thin, supple material.



Pro-ophta® Junior

Skin coloured nonwoven carrier: 100% polyester; occlusive air-permeable centre: fabric 100% cotton; hypoallergenic, porously coated synthetic rubber adhesive.



Pro-ophta® Junior Mini

6.5 x 5.4 cm	REF	Shipping Units
in a plastic bag	16 841	5 pcs. / 2000 pcs.
in a cardboard box	34 231	20 pcs. / 800 pcs.
in a cardboard box	16 846	50 pcs. / 1300 pcs.
in a cardboard box	16 843	100 pcs. / 2400 pcs.

Pro-ophta® Junior Maxi

7.0 x 5.9 cm	REF	Shipping Units
in a plastic bag	16 840	5 pcs. / 2000 pcs.
in a cardboard box	16 844	20 pcs. / 800 pcs.
in a cardboard box	16 845	50 pcs. / 1300 pcs.
in a cardboard box	16 842	100 pcs. / 2400 pcs.

Ready-for-use dressings made from widthwise elastic fabric; for the treatment of injuries to the external eye and for protection after eye operations.

Pro-ophta eye dressings mould themselves to the contours of the face. The hypoallergenic adhesive ensures reliable retention even in long-term treatment. The dressings are nevertheless easy to remove without leaving any residues.

Pro-ophta® Eye Dressings

Widthwise elastic fabric, skin coloured, 70% viscose, 30% polyamide, **transparent disk:** clear polycarbonate, **occlusive cover:** black polycarbonate, **mesh/gauze pad:** gauze 100% cotton; mesh 100% viscose; hypoallergenic polyacrylate adhesive.



Pro-ophta® Eye Dressing S

Elastic dressing with transparent, shaped disk. Mainly used following ptosis operations and lagophthalmus as well as facial palsy. Available in two sizes to suit various groups of patients.

Pro-ophta® Eye Dressing D

Elastic dressing with occlusive, shaped disc. To protect the eye following operations, in particular after ptosis operations, in cases of facial palsy as well as lagophthalmus, also to cover the eye during the night.

Pro-ophta® Eye Dressing K

Elastic dressing with incorporated cotton gauze pad; for the treatment of injuries to the external eye, also to immobilize the eyelid.

Pro-ophta® Eye Dressing S			
individually sealed		REF	Shipping Units
non-sterile	small	34 227	1 pc. / 200 pcs.
		34 228	50 pcs. / 600 pcs.
	large	34 229	1 pc. / 120 pcs.
		34 230	50 pcs. / 600 pcs.
Pro-ophta® Eye Dressing D			
individually sealed			
non-sterile		34 225	1 pc. / 200 pcs.
		34 226	50 pcs. / 600 pcs.
Pro-ophta® Eye Dressing K			
individually sealed			
non-sterile		34 104	50 pcs. / 600 pcs.