

Pro-ophta® Protective Eye Padding Ring

Extremely soft special padding ring for postoperative eye care which maintains the natural eyelid function, especially after cataract operations, eyelid surgery and injuries to the external eye.

Pro-ophta protective padding rings are worn under perforated eye shields. The soft, smooth padding rings are lint-free and absorbent. They provide cushioning and allow plenty of air to reach the eye.

The protective eye padding rings with perforated eye shields are popular with patients because they permit early testing of vision – a psychologically important factor for patient compliance.



Pro-ophta® Protective Eye Padding Ring

Soft nonwoven, 67% viscose, 33% polyester, in 16 layers joined by quilted seams and punched to the required ring shape with open centre.



The perforated plastic eye shield is not part of the product.

Pro-ophta® Protective Eye Padding Ring		REF	Shipping Units
sterile individually sealed	7 x 8.5 cm	13 014	50 pcs. / 500 pcs.
non-sterile in a cardboard box	7 x 8.5 cm	13 015	50 pcs. / 300 pcs.