

# Pro-ophta<sup>®</sup> Sticks

## Pro-ophta<sup>®</sup> Sticks

Tissue-preserving instrument made from soft, binder-free nonwoven, rolled into sticks and obliquely cut at both ends. The sticks are used for tissue preparation and swabbing in ophthalmology.

Pro-ophta sticks have a good wicking action but retain the strength required for swabbing and preparing.

Highly absorbent nonwoven, white, rolled into sticks, strengthened externally with a fine coverstock.



### Pro-ophta<sup>®</sup> Sticks

ø 5 mm, 66 mm long

REF

Shipping Units

**sterile**, 20 sticks packed in a peel pouch  
**non-sterile**

16 516  
16 515

400 pcs. / 1600 pcs.  
500 pcs. / 25000 pcs.