



# Raucocel® Ear Packings

Once expanded, Raucocel ear packings conform exactly and guarantee, through their elasticity, optimum medical supply to

the auditory canal lumen, even while the swelling decreases. The material's high capillarity prevents the inserted solution

from sinking down.



## Ear Packing without drawstring


mainly for the treatment of inflammations of the external auditory canal.



## Ear Packing with drawstring

for the postoperative treatment, e.g. after tympanoplasty.

The packing keeps the auditory canal open during the healing process and keeps skin transplants in place. With drawstring.



Ear Packing sterile	REF	Shipping Units
9 x 15 mm, without drawstring	10 171	50 pcs.
9 x 15 mm, with drawstring	10 172	10 pcs.
12 x 15 mm, with drawstring	10 169	10 pcs.