



Raucocel® Epistaxis Packings

Epistaxis Packings work with the principle of autocoagulation and thereby support the physiological clotting.

A specific benefit is the time-saving and atraumatic application of the packing. **Before removal the packing must be thoroughly hydrated;** this virtually prevents the risk of after-bleeding.


Epistaxis Packing 100

for the control of posterior nasal bleeding and for use in postoperative care, e.g. after conchotomies.




Epistaxis Packing 55

for the control of medial nasal bleeding.




Epistaxis Packing 100 sterile		
	REF	Shipping Units
100 mm, with drawstring	10 143	10 pcs. / 240 pcs.



10 143

Epistaxis Packing 55 sterile		
	REF	Shipping Units
55mm, with drawstring	10 144	20 pcs. / 240 pcs.



10 144

Raucocel® Accessories		
	REF	Shipping Units
Instrument Wipe, sterile 76 x 76 mm	10 313	20 pcs. / 480 pcs.
Headband	11 971	10 pcs. / 240 pcs.
Diamond Cleaning Block, sterile	11 483	10 pcs. / 120 pcs.