


Raucocel® Nasal Packings


Straight Nasal Packing

recommended for use after rhinoplasties. The most simple version of postoperative nasal packings has a straight shape and no airway tube. The version with airway tube provides an improved equalisation of pressure and eases swallowing. At the same time the airway tube reinforces the packing's stability. With drawstring.




Anatomically shaped Nasal Packing

guarantees optimum support after septo-/rhinoplasties, for small noses. The version with airway tube provides better equalisation of pressure and eases swallowing. At the same time the airway tube reinforces the packing's stability. With drawstring.




Slim Nasal Packing

recommended after septum- and septo-/rhinoplasties, for small noses. The anatomically shaped Slim Nasal Packing has only 2/3 of the width of the anatomically shaped packing and is ideally suited for smaller noses. With drawstring.




Sinus Nasal Packing "K"

for the postoperative treatment after endoscopic sinus surgery, especially after ethmoidectomy. Its special shape makes it particularly suitable for tamponing the ethmoidal sinus to avoid deformities. With drawstring.




Straight Nasal Packing sterile

	REF	Shipping Units
45 mm, with airway tube, with drawstring	11 004	10 pcs. / 240 pcs.
80 mm, with airway tube, with drawstring	11 011	10 pcs. / 240 pcs.
45 mm, with drawstring	11 002	10 pcs. / 240 pcs.
80 mm, with drawstring	11 010	10 pcs. / 240 pcs.

**Anatomically shaped Nasal Packing sterile**

	REF	Shipping Units
55 mm, with airway tube, with drawstring	11 014	10 pcs. / 240 pcs.
80 mm, with airway tube, with drawstring	11 015	10 pcs. / 240 pcs.

**Sinus Nasal Packing "K" sterile**

	REF	Shipping Units
9 mm, with drawstring	11 937	20 pcs. / 480 pcs.
12 mm, with drawstring	11 935	20 pcs. / 480 pcs.

