


Raucodrape® Ophthalmic Drapes

Raucodrape® Surgical Incision Drape

Raucodrape® Ophthalmic Drape 150 x 150 cm

Patient drape with central aperture with integrated incision film and fluid collection pouch on one side. For ophthalmological operations, e.g. cataract surgery.



150 x 150 cm

Raucodrape® Surgical Incision Drape

Elastic and supple polyurethane film, coated with polyacrylate adhesive, hypoallergenic, prevents the migration of micro-organisms, but is permeable to moisture vapour, transparent and non-reflective. Raucodrape film incision drape is used for covering the operating site.

Polyurethane film, elastic and transparent, coated with polyacrylate adhesive, sterilised by radiation.

Raucodrape® Ophthalmic Drapes sterile			
	REF	Shipping Units	
150 x 150 cm, with aperture 7 x 10 cm, integrated incision film and fluid collection pouch on one side	33 056	17 pcs. / 34 pcs.	
additionally for ophthalmological operations: Armrest Cover, in pairs, 35 x 73 cm	33 123	46 pcs. / 184 pcs.	
Raucodrape® Surgical Incision Drape sterile, individually sealed			
adhesive surface	12 x 17 cm	25 440	10 pcs. / 80 pcs.
	15 x 20 cm	25 441	10 pcs. / 60 pcs.
	30 x 20 cm	25 442	10 pcs. / 80 pcs.
	45 x 20 cm	25 443	10 pcs. / 80 pcs.