

First Aid Products

First Aid Cabinet

Plastic, approx. 33 x 45 x 16 cm, with adjustable plastic shelves and storage space inside the door; one door, white, without contents.

Small Refill

DIN 13157-C

Contents¹⁾

- 1 Porofix First-aid dressing
DIN 13019-A, 2,5 cm x 5 m
- 8 Curaplast Sensitive adhesive wound dressings DIN 13019-E, 6 x 10 cm
- 5 Finger tip dressings
- 5 Curaplast Sensitive finger dressings
DIN 13 019-E, 2 x 18 cm
- 10 Curaplast Sensitive wound dressings
1.9 x 7.2 cm
- 3 Dressing kits DIN 13151-M, sterile
- 2 Dressing kits DIN 13151-G, sterile
- 1 Metalline dressing cloth
DIN 13152-A, 60 x 80 cm
- 1 Metalline dressing cloth
DIN 13152-BR, 40 x 60 cm
- 6 Gazin wound dressings
10 x 10 cm, sterile, in pairs
- 2 Eye dressings 5.5 x 7.5 cm, sterile
- 1 Aluminized film, 160 x 220 cm
- 3 Conforming bandages DIN 61634,
8 cm
- 3 Conforming bandages DIN 61634,
6 cm
- 1 tg fix tubular net bandage, size E,
4 m stretched
- 1 Triangular sheet DIN 13168-D
- 1 Pair of scissors DIN 58279
- 10 Nonwoven cloths, 20 x 30 cm
- 2 PE-bags, sealable
- 4 Disposable gloves
- 1 Brief information on immediate
measures to be taken at the site of
an accident
List of contents

Large Refill

DIN 13169-E

Contents¹⁾

- 1 Porofix First-aid dressing
DIN 13019-A, 2,5 cm x 5 m
- 16 Curaplast Sensitive adhesive wound dressings DIN 13019-E, 6 x 10 cm
- 10 Finger tip dressings
- 10 Curaplast Sensitive finger dressings
DIN 13 019-E, 2 x 18 cm
- 20 Curaplast Sensitive wound dressings
1.9 x 7.2 cm
- 6 Dressing kits DIN 13151-M, sterile
- 4 Dressing kits DIN 13151-G, sterile
- 2 Metalline dressing cloths
DIN 13152-A, 60 x 80 cm
- 2 Metalline dressing cloths
DIN 13152-BR, 40 x 60 cm
- 12 Gazin wound dressings
10 x 10 cm, sterile, in pairs
- 4 Eye dressings 5.5 x 7.5 cm, sterile
- 2 Aluminized films, 160 x 220 cm
- 6 Conforming bandages DIN 61634,
8 cm
- 6 Conforming bandages DIN 61634,
6 cm
- 1 tg fix tubular net bandage, size E,
4 m stretched
- 1 tg fix tubular net bandage, size F,
4 m stretched
- 2 Triangular sheets DIN 13168-D
- 1 Pair of scissors DIN 58279
- 20 Nonwoven cloths, 20 x 30 cm
- 4 PE-bags, sealable
- 8 Disposable gloves
- 1 Brief information on immediate
measures to be taken at the site of an
accident
List of contents

¹⁾Subject to alterations

Triangular bandage		REF	Shipping Units
DIN 13168-D, sealed in a foil pouch, unbleached white, 96 x 96 x 136 cm		14 970	1 pc. / 250 pcs.
Refills			
small		43 103	1 pc.
large		43 104	1 pc.
Dressing Kit, sterile			
DIN 13151, sealed in a foil pouch, impermeable to air and water	small	14 900	1 pc. / 650 pcs.
	medium	14 901	1 pc. / 500 pcs.
	large	14 902	1 pc. / 500 pcs.