



Suprasorb[®] CNP P2

For the mobile patient

The L&R Controlled Negative Pressure therapy offers the right solution for the specific needs of your mobile patients.

Comfortable wear and ergonomic design are prerequisites for patient compliance and thereby decisive factors leading to your healing success.



The Suprasorb® CNP P2 negative pressure pump is highly suitable for mobile patients who can leave their bed by themselves, or for outpatient care.

The mobility gained improves the quality of life and patient compliance. Early mobilisation and the resulting transfer to outpatient care reduce the cost of treatment.

Convincing benefits:

- low weight and size
- capacity 250 ml
- easy, mobile handling through operation in the carrier bag provided
- comfortable to wear through ergonomic design
- carrier bag with adjustable strap, can be worn over the shoulder or on the hip
- internal bacteria filter with integrated overflow and odour control
- exceptionally quiet, with low vibration
- closed exudate collector with filling level easy to read from outside
- operation possible with the unit lying down
- all disposable parts are easy to exchange

User-friendly

The inspection window in the bag allows permanent monitoring of the therapy and uncomplicated operation, even when the bag is closed



The Suprasorb® CNP P2 pump is anatomically designed and therefore comfortable to wear – both over the shoulder or on the hip.

In combination with the low weight and the extremely quiet operation, the patient can participate in everyday life almost without limitations.

